



ELISA KITS CATALOG

ENZYME-LINKED IMMUNOSORBENT ASSAY KIT

Your Good Partner In Biology Research




Feiyue Biotechnology Co., LTD
Your Good Partner In Biology Research
Contact us
Tel: +86 027 8700 2654
E-mail: info@feiyuebio.com

Address:
Building B3, No.666, Gaoxin Avenue, East Lake New Technology Development Zone, Wuhan, Hubei Province, China

FEIYUEBIO

ELISA KITS

Research Area:
Our research area is Cytokines, Cancer, Cardiovascular,
Immunology, Tumor marker, Growth Factors/Hormones,
TORCH..etc.




About Us

FEIYUEBIO As your biological research partner, we offer 5000+ ELISA test kits, primarily in sandwich ELISA format, covering species such as human, mouse, rat, and bovine...

As well as an accurate research tool for the detection and quantification of targets of interest in cultures and samples. ELISA kits are used in a variety of important research fields, including cancer, cardiovascular, cell biology, epigenetics, metabolism, immunology, neuroscience, signal transduction, and so on.

FEIYUEBIO maintains the strictest quality control system to ensure that every ELISA kit is free of flaws, allowing you to obtain accurate data in your experiments. We offer a wide-ranging portfolio of human and animal ELISA kits that are ultra-sensitive and can be performed in as little as 90 minutes.



IMMUNOLOGY

FeiyueBio has a large number of ELISA Kits related to immunology research, including:

1. High sensitivity, can detect extremely low concentration of cytokines
2. Covering a wide range, involving the main research directions of immunology, The sample volume of mice, rats and other species is as small as 10-50ul,

METABOLISM

FeiyueBio Metabolism-related ELISA kit has:


1. Covering a wide range of species, including chickens, pigs, cattle, sheep, bovine, hamsters and other rare species...
2. Focus on the classic fat metabolism-related cytokines related to obesity, such as: Adiponectin/ADPN, Leptin, Resistin, ghrelin, etc.
3. Strict quality control to ensure the stability of the kit.

Research Areas

CANCER

FeiyueBio Cancer related ELISA kit has:


1. Covering a wide range of species, in addition to common species such as humans, mice, rats, and rabbits, there are also rare species such as chickens, dogs, and hamsters...
2. Focus on some newer but promising tumor markers, such as: pepsinogen I/II related to gastric cancer.



SIGNAL TRANSDUCTION

The related ELISA kit of FeiyueBio Signal Transduction has:

1. The linearity is good, and the R square value of each kit standard curve reaches more than 29.
2. The detection range is wide, and the dilution factor of most samples is small Involves a wide range, including JAK-STAT pathway, NF-κB pathway, Toll-like receptor pathway, TGF-β pathway, p53 pathway, etc.



CARDIOVASCULAR

FeiyueBio Cardiovascular related ELISA kit has:


1. Strong specificity, no cross-reaction with other analogues.
2. Strict quality inspection ensures the minimum intra-assay and inter-assay differences of the kit to the greatest extent
3. Focus on potential new indicators, such as: the best indicator of heart failure--soluble growth-stimulating expression gene 2 protein (ST2).



HEMATOPATHY

FeiyueBio Hematopathy related ELISA kit has:

1. Covering a wide range of species, in addition to common species, it also involves rare species such as chickens, pigs, dogs, hamsters, etc.
2. Strict quality control to ensure the stability of the system and good linearity of the standard curve.
3. Good specificity and minimal cross-reaction.



Immunology Series Products

Cata NO.	Product Name	Detection range	Sensitivity
FY-EU1002	PGE1(Prostaglandin E1) ELISA Kit	31.25-2000 pg/mL	9.55 pg/mL
FY-EU1023	FK506(Tacrolimus) ELISA Kit	0.16-10 ng/mL	0.053ng/mL
FY-EM1157	Mouse ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	62.5-4000ng/mL	37.5 ng/mL
FY-EU1211	Testo(Testosterone) ELISA Kit	0.16-10ng/mL	0.09 ng/mL
FY-EU1219	IS(serum indoxyl sulfate) ELISA Kit	31.25-2000ng/mL	18.75 ng/mL
FY-EU1222	8-OHdG(8-Hydroxydeoxyguanosine) ELISA Kit	1.56~100 ng/mL	0.94 ng/mL
FY-EH1266	Human DKK1(Dickkopf Related Protein 1) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-ER1272	Rat EGF(Epidermal Growth Factor) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1281	Human ET-1(Endothelin 1) ELISA Kit	1.25-80pg/mL	0.75 pg/mL
FY-EM1324	Mouse FST(Follistatin) ELISA Kit	0.31-20ng/mL	0.19 ng/mL
FY-EH1405	Human TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	7.81-500pg/mL	3.91 pg/mL
FY-EH1414	Human Bcl-2(B-cell Leukemia/Lymphoma 2) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM1445	Mouse IL-1β(Interleukin 1 Beta) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM1463	Mouse TE(telomerase) ELISA Kit	6.25-400ng/mL	3.75 ng/mL
FY-EM1514	Mouse NCAD(Neural Cadherin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1570	Human ADA(Adenosine Deaminase) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1649	Human ALDM(Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EU1675	LPS(Lipopolysaccharides) ELISA Kit	1.57-100 ng/mL	0.68 ng/mL
FY-EH1801	Human LCN1(Lipocalin 1) ELISA Kit	0.79-50ng/mL	0.47 ng/mL
FY-EH1864	Human ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-ER1874	Rat AST(Aspartate Aminotransferase) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1877	Human ATM(Ataxia Telangiectasia Mutated) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1944	Human MT1E(Metallothionein 1E) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EM1989	Mouse CALU(Calumenin) ELISA Kit	7.82-500pg/mL	4.69 pg/mL
FY-EH2005	Human cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EM2068	Mouse CFL1(Cofilin 1, Non-Muscle) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH2095	Human MUC6(Mucin 6) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH2096	Human MYBPC3(Myosin Binding Protein C, Cardiac) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-ER2170	Rat CFL2(Cofilin 2, Muscle) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM2200	Mouse C3a(Complement Component 3a) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM2260	Mouse GATA4(GATA Binding Protein 4) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EM2437	Mouse GAL12(Galectin 12) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EM2440	Mouse GAL2(Galectin 2) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM2444	Mouse GAL4(Galectin 4) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM2486	Mouse FOXP3(Forkhead Box Protein P3) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EM2539	Mouse ELN(Elastin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH2653	Human SEMA3A(Semaphorin 3A) ELISA Kit	0.16-10ng/mL	0.10 ng/mL
FY-EM2731	Mouse GR(glucocorticoid receptor) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH2938	Human PCT(Procalcitonin) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM2947	Mouse CNTN1(Contactin 1) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM3058	Mouse VEGFR-2/KDR ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-ER3139	Rat CX43(Connexin 43) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH3159	Human Integrin beta-1(ITGB1) ELISA kit	62.5-4000 pg/mL	15.6 pg/mL
FY-EH3170	Human GPC4(Glypican 4) ELISA Kit	0.16-10ng/mL	0.10 ng/mL
FY-EH3172	Human TRUB1(Probable tRNA Pseudouridine Synthase 1) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL

Cancer Series Products

Cata NO.	Product Name	Detection Range	Sensitivity
FY-EU1013	TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1014	AcCoA(Acetyl CoA) ELISA Kit	0.32-20 nmol/mL	0.119 nmol/mL
FY-EM1044	Mouse BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EU1077	TMAO(Trimethylamine-N-oxide) ELISA Kit	15.63-1000 ng/mL	4.93 ng/mL
FY-EU1129	Cor(Cortisol) ELISA Kit	0.31-20ng/mL	0.19 ng/mL
FY-EU1141	Pg(Progesterone) ELISA Kit	78.13-5000pg/mL	46.88 pg/mL
FY-EU1211	Testo(Testosterone) ELISA Kit	0.16-10ng/mL	0.09 ng/mL
FY-EU1213	E2(Estradiol) ELISA Kit	23.44-1500 pg/mL	14.06 pg/mL
FY-EU1222	8-OHdG(8-Hydroxydeoxyguanosine) ELISA Kit	1.56~100 ng/mL	0.94 ng/mL
FY-EM1267	Mouse DKK1(Dickkopf Related Protein 1) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EM1324	Mouse FST(Follistatin) ELISA Kit	0.31-20ng/mL	0.19 ng/mL
FY-EZe1339	Zebrafish VTG(Vitellogenin) ELISA Kit	31.25-2000 pg/mL	9.34 pg/mL
FY-EH1405	Human TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	7.81-500pg/mL	3.91 pg/mL
FY-EH1436	Human PRL(Prolactin) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1444	Human IL-1β(Interleukin 1 Beta) ELISA Kit	7.82-500pg/mL	4.69 pg/mL
FY-EM1445	Mouse IL-1β(Interleukin 1 Beta) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-ER1449	Rat IL-12(Interleukin 12) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM1480	Mouse COL4(Collagen Type IV) ELISA Kit	1.57-100ng/mL	0.94 ng/mL
FY-EM1491	Mouse PⅢNP(N-Terminal Procollagen Ⅲ Propeptide) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1498	Human AGGF1 ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM1514	Mouse NCAD(Neural Cadherin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1517	Human PDGFRL ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1649	Human ALDM(Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1765	Human ALOXE3(Arachidonate Lipooxygenase 3) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1775	Human IGF1R(Insulin Like Growth Factor 1 Receptor) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1783	Human LTβR(Lymphotoxin Beta Receptor) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EM1829	Mouse AIF(Apoptosis Inducing Factor) ELISA Kit	7.82-500pg/mL	4.69 pg/mL
FY-EM1834	Mouse ASK1(Apoptosis Signal Regulating Kinase 1) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-ER1902	Rat BAFFR(B-Cell Activation Factor Receptor) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-ER1912	Rat Bcl-10(B-cell Lymphoma/Leukemia 10) ELIS Kit ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1969	Human CTR(Calcitonin Receptor) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM1983	Mouse CRT(Calreticulin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH2005	Human cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EM2048	Mouse CTSK(Cathepsin K) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH2060	Human CCAAT/Enhancer Binding Protein, Epsilon ELISA Kit	3.13-200ng/mL	1.88 ng/mL
FY-EH2143	Human NOS3/eNOS ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EM2161	Mouse F5(Coagulation Factor V) ELISA Kit	3.13-200ng/mL	1.88 ng/mL
FY-ER2526	Rat TRH(Thyrotropin Releasing Hormone) ELISA Kit	7.82-500pg/mL	4.69 pg/mL
FY-EM2528	Mouse GHRH(Growth Hormone Releasing Hormone) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM2542	Mouse ELA4(Elastase 4) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH2653	Human SEMA3A(Semaphorin 3A) ELISA Kit	0.16-10ng/mL	0.10 ng/mL
FY-ER2720	Rat OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM2752	Mouse ON(Osteonectin) ELISA Kit	3.13-200ng/mL	1.88 ng/mL
FY-EM2765	Mouse Peroxisome Proliferator Activated Receptor Alpha ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM2768	Mouse PPAR-γ ELISA Kit	0.16-10ng/mL	0.1 ng/mL

Metabolism Series Products

Cata NO.	Product Name	Detection Range	Sensitivity
FY-EU1004	F-E3(Free Estriol) ELISA Kit	0.63-40 ng/mL	0.218 ng/mL
FY-EU1005	F-TESTO(Free Testosterone) ELISA Kit	31.25-2000 pg/mL	18.75 pg/mL
FY-EU1006	VA(Vitamin A) ELISA KitKit	15.63-1000 ng/mL	9.38 ng/mL
FY-EU1007	4-HNE(4-Hydroxynonenal) ELISA Kit	0.63-40 ng/mL	0.38 ng/mL
FY-EU1009	SP(Substance P) ELISA KitKit	78.13-5000 pg/mL	46.88 pg/mL
FY-EU1013	TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1015	25-OH-D(25 Hydroxy Vitamin D) ELISA Kit	3.13-200 ng/mL	0.91 ng/mL
FY-EU1016	FSH(Follicle Stimulating Hormone) ELISA Kit	1.57-100mIU/mL	0.67 mIU/mL
FY-EU1021	Cit(Citrulline) ELISA Kit	1.57-100 ng/mL	0.51 ng/mL
FY-EU1024	CTN(Creatinine) ELISA Kit	3.13-200 μg/mL	0.89 μg/mL
FY-EU1025	11-KT(11-keto Testosterone) ELISA Kit	156.25-10000 pg/mL	15.29 pg/mL
FY-EU1026	ACTb(Actin Beta) ELISA Kit	0.16-10 ng/mL	0.057 ng/mL
FY-EM1030	Mouse BAFF/CD257(B-cell Activating Factor) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EU1036	DOPAC (3,4-dihydroxyphenyl acetic acid) Elisa Kit.	1.57-100 ng/mL	0.59 ng/mL
FY-EU1037	DHEA-S(Dehydroepiandrosterone sulfate) ELISA Kit	6.25-400 ng/mL	4.25 ng/mL
FY-EU1039	CH(Cholesterol) ELISA Kit	7.82-500 μg/mL	2.48 μg/mL
FY-EU1040	SMD(Spermidine) ELISA Kit	3.13-200 ng/mL	0.87 ng/mL
FY-EU1041	VB7(Vitamin B7) ELISA Kit	312.5-20000 pg/mL	90.2 pg/mL
FY-EH1043	Human BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EU1046	CS(Chondroitin Sulfate) ELISA Kit	156.25-10000 pg/mL	15.17 pg/mL
FY-EU1053	VB2(Vitamin B2) ELISA Kit	9.38-600 ng/mL	2.55 ng/mL
FY-EH1055	Human BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EM1056	Mouse BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-ER1057	Rat BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EU1064	VMA(Vanillylmandelic Acid) ELISA Kit	156.25-10000 ng/mL	43.7 ng/mL
FY-EU1067	CP(Cloprostenol) ELISA Kit	156.25-10000 pg/mL	44.3 pg/mL
FY-EH1068	Human BMP-7(Bone Morphogenetic Protein 7) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1076	ALCAR(Acetylcarnitine) ELISA Kit	78.13-5000 ng/mL	25.2 ng/mL
FY-EU1085	PGH2(Prostaglandin H2) ELISA Kit	3.13-200 pg/mL	0.97 pg/mL
FY-EU1086	TG(Triglyceride) ELISA Kit	15.63-1000 μg/mL	1.8 μg/mL
FY-EU1088	LNG(Levonorgestrel) ELISA Kit	1.57-100 ng/mL	0.53 ng/mL
FY-EU1089	12-HETE(12-Hydroxyeicosatetraenoic Acid) ELISA Kit	15.63-1000 ng/mL	5.32 ng/mL
FY-EU1091	PS(Phosphatidylserine) ELISA Kit	18.75-1200 ng/mL	5.27 ng/mL
FY-EU1092	IP3(Inositol Triphosphate) ELISA Kit	156.25-10000 pg/mL	49.3 pg/mL
FY-EH1093	Human CA9(Carbonic Anhydrase IX) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EU1096	DAG(Diacylglycerol) ELISA Kit	1562.5-100000 pg/mL	141 pg/mL
FY-EU1102	VK1(Vitamin K1) ELISA Kit	156.25-10000 pg/mL	49.3 pg/mL
FY-EU1104	VC(Vitamin C) ELISA Kit	625-40000 ng/mL	192.1 ng/mL
FY-EU1112	PGI2(Prostacyclin) ELISA Kit	3.13-200 pg/mL	0.93 pg/mL
FY-EU1114	VD3(Vitamin D3) ELISA Kit	6.25-400 ng/mL	2.14 ng/mL
FY-EU1116	DHVD3(1,25-Dihydroxyvitamin D3) ELISA Kit	156.25-10000 pg/mL	49.5 pg/mL
FY-EH1117	Human MIP-1α ELISA Kit	23.44-1500pg/mL	14.06 pg/mL
FY-EU1120	PGD2(Prostaglandin D2) ELISA Kit	15.63-1000 pg/mL	4.34 pg/mL
FY-EU1121	Car(Carnosine) ELISA Kit	14.07-900 ng/mL	3.53 ng/mL
FY-EU1122	AED(Androstenediol) ELISA Kit	156.25-10000 pg/mL	53.4 pg/mL

Signal Transduction Series Products

Cata NO.	Product Name	Detection Range	Sensitivity
FY-EU1007	4-HNE(4-Hydroxynonenal) ELISA Kit	0.63-40 ng/mL	0.38 ng/mL
FY-EU1008	GnRH(GonadotropinReleasing Hormone) ELISA Kit	15.63-1000 pg/mL	9.38 pg/mL
FY-EU1010	MEK(Methionine-Enkephalin) ELISA Kit	125.00-8000 pg/mL	75.00 pg/mL
FY-EU1013	TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1014	AcCoA(Acetyl CoA) ELISA Kit	0.32-20 nmol/mL	0.119 nmol/mL
FY-EM1018	Mouse ANG2(Angiopoietin 2) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1024	CTN(Creatinine) ELISA Kit	3.13-200 µg/mL	0.89 µg/mL
FY-EU1026	ACTb(Actin Beta) ELISA Kit	0.16-10 ng/mL	0.057 ng/mL
FY-EH1029	Human BAFF/CD257(B-Cell Activating Factor) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EU1032	MGO(Methylglyoxal) ELISA Kit	1.57-100 ng/mL	0.51 ng/mL
FY-EU1033	HS(Heparan Sulfate) ELISA Kit	0.31-20 ng/mL	0.19 ng/mL
FY-EM1044	Mouse BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EU1061	PIP3(Phosphatidylinositol Trisphosphate) ELISA Kit	7.82-500 ng/mL	2.83 ng/mL
FY-EU1062	GS(Guanosine) ELISA Kit	312.5-20000 pg/mL	82.7 pg/mL
FY-EU1066	PI(Phosphatidylinositol) ELISA Kit	156.25-10000 pg/mL	48.6 pg/mL
FY-ER1070	Rat BMP-7(Bone Morphogenetic Protein 7) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EU1072	Tau(Taurine) ELISA Kit	15.63-1000 ng/mL	4.39 ng/mL
FY-EU1078	ADP(Adenosine Diphosphate) ELISA Kit	23.44-1500 ng/mL	7.17 ng/mL
FY-EU1079	BA(Butyric Acid) ELISA Kit	156.25-10000 pg/mL	51.5 pg/mL
FY-EU1080	TCA(Taurocholic Acid) ELISA Kit	4.69-300 ng/mL	1.41 ng/mL
FY-EU1084	SPH(Sphingomyelin) ELISA Kit	6.25-400 µg/mL	1.85 µg/mL
FY-EU1089	12-HETE(12-Hydroxyeicosatetraenoic Acid) ELISA Kit	15.63-1000 ng/mL	5.32 ng/mL
FY-EU1092	IP3(Inositol Triphosphate) ELISA Kit	156.25-10000 pg/mL	49.3 pg/mL
FY-EU1096	DAG(Diacylglycerol) ELISA Kit	1562.5-100000 pg/mL	141 pg/mL
FY-EU1098	cGMP(Cyclic Guanosine Monophosphate) ELISA Kit	0.78-50 pmol/mL	0.47 pmol/mL
FY-EU1109	cAMP(Cyclic Adenosine Monophosphate) ELISA Kit	312.5-20000 pg/mL	92.5 pg/mL
FY-EU1129	Cor(Cortisol) ELISA Kit	0.31-20ng/mL	0.19 ng/mL
FY-EU1133	Nicotinic Acid Adenine Dinucleotide Phosphate ELISA Kit	156.25-10000 ng/mL	48.7 ng/mL
FY-EU1141	Pg(Progesterone) ELISA Kit	78.13-5000pg/mL	46.88 pg/mL
FY-EU1155	ACH(Acetylcholine) ELISA Kit	15.63-1000 pg/mL	5.31 pg/mL
FY-EU1161	NE(Noradrenaline) ELISA Kit	78.13-5000 pg/mL	26.1 pg/mL
FY-EU1164	DA(Dopamine) ELISA Kit	15.63-1000 pg/mL	4.71 pg/mL
FY-EU1174	HA(Histamine) ELISA Kit	1.57-100 ng/mL	0.51 ng/mL
FY-EU1184	IAA(Indole 3 Acetic Acid) ELISA Kit	3.13-200 ng/mL	0.89 ng/mL
FY-EU1185	LXA4(Lipoxin A4) ELISA Kit	0.78-50 ng/mL	0.47 ng/mL
FY-EU1187	TXA2(Thromboxane A2) ELISA Kit	15.63-1000 pg/mL	4.39 pg/mL
FY-EU1191	AA(Arachidonic Acid) ELISA Kit	3.13-200 µg/mL	0.94 µg/mL
FY-EM1194	Mouse MIPF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-ER1195	Rat MIPF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1196	bHB(Beta-Hydroxybutyric Acid) ELISA Kit	1.57-100 µg/mL	0.48 µg/mL
FY-EU1202	5-HT(5-Hydroxytryptamine) ELISA Kit	14.07-900 ng/mL	3.97 ng/mL
FY-EU1204	ft4(Free Thyroxine) ELISA Kit	1.41-90 pg/mL	0.48 pg/mL
FY-EU1205	ft3(Free Triiodothyronine) ELISA Kit	1.56-100 pg/mL	0.94 pg/mL
FY-EU1209	gABA(Gamma-Aminobutyric Acid) ELISA Kit	31.25-2000 pg/mL	9.27 pg/mL
FY-EU1211	Testo(Testosterone) ELISA Kit	0.16-10ng/mL	0.09 ng/mL

Cardiovascular Series Products

Cata NO.	Product Name	Detection Range	Sensitivity
FY-EU1013	TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1017	Human ANG2(Angiopoietin 2) ELISA Kit	46.88-3000pg/mL	28.13 pg/mL
FY-EU1064	VMA(Vanillylmandelic Acid) ELISA Kit	156.25-10000 ng/mL	43.7 ng/mL
FY-EU1075	LA(Lactic Acid) ELISA Kit	9.38-600 µg/mL	2.92 µg/mL
FY-EU1077	TMAO(Trimethylamine-N-oxide) ELISA Kit	15.63-1000 ng/mL	4.93 ng/mL
FY-EH1081	Human E-Cad(E-Cadherin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EU1086	TG(Triglyceride) ELISA Kit	15.63-1000 µg/mL	1.8 µg/mL
FY-EU1092	IP3(Inositol Triphosphate) ELISA Kit	156.25-10000 pg/mL	49.3 pg/mL
FY-EU1112	PGI2(Prostacyclin) ELISA Kit	3.13-200 pg/mL	0.93 pg/mL
FY-EU1124	PAF(Platelet Activating Factor) ELISA Kit	46.88-3000 pg/mL	14.4 pg/mL
FY-EU1129	Cor(Cortisol) ELISA Kit	0.31-20ng/mL	0.19 ng/mL
FY-EH1130	Human MIP-1β(Macrophage Inflammatory Protein 1 Beta) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1141	Pg(Progesterone) ELISA Kit	78.13-5000pg/mL	46.88 pg/mL
FY-EU1151	CML(Carboxymethyl Lysine) ELISA Kit	78.13-5000 ng/mL	25.9 ng/mL
FY-EU1164	DA(Dopamine) ELISA Kit	15.63-1000 pg/mL	4.71 pg/mL
FY-EH1169	Human TARC(Thymus Activation Regulated Chemokine) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1177	DDAVP(1-Desamino 8D Arginine Vasopressin) ELISA Kit	156.25-10000 pg/mL	39.5 pg/mL
FY-EH1181	Human SLC(Secondary Lymphoid Tissue Chemokine) ELISA Kit	78.13-5000pg/mL	46.88 pg/mL
FY-EU1186	PGF2a(Prostaglandin F2 Alpha) ELISA Kit	12.5-800 pg/mL	3.88 pg/mL
FY-EU1189	ADMA(Asymmetrical Dimethylarginine) ELISA Kit	15.63-1000 ng/mL	4.51 ng/mL
FY-EH1193	Human MIPF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	78.13-5000pg/mL	46.88 pg/mL
FY-EH1206	Human MIP4α(Macrophage Inflammatory Protein 4 Alpha) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1207	AGE(Advanced Glycation End Product) ELISA Kit	125-8000 ng/mL	38.2 ng/mL
FY-EU1209	gABA(Gamma-Aminobutyric Acid) ELISA Kit	31.25-2000 pg/mL	9.27 pg/mL
FY-EU1220	PROA (Protein A) ELISA Kit	46.88-3000 pg/mL	18.6 pg/mL
FY-EU1222	8-OHdG(8-Hydroxydeoxyguanosine) ELISA Kit	1.56~100 ng/mL	0.94 ng/mL
FY-EM1243	Mouse GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EGP1253	Guinea Porcine ADMA(Asymmetrical Dimethylarginine) ELISA Kit	15.63-1000 ng/mL	4.51 ng/mL
FY-ER1258	Rat BLC(B-Lymphocyte Chemoattractant) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM1264	Mouse Cyt-C(Cytochrome C) ELISA Kit	0.4-25ng/mL	0.24 ng/mL
FY-EZe1273	Zebrafish DA(Dopamine) ELISA Kit	15.63-1000 pg/mL	4.71 pg/mL
FY-ER1312	Rat FGF21(Fibroblast Growth Factor 21) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-ECH1315	Chicken DA(Dopamine) ELISA Kit	15.63-1000 pg/mL	4.71 pg/mL
FY-EP1322	Porcine IP3(Inositol Triphosphate) ELISA Kit	156.25-10000 pg/mL	49.7 pg/mL
FY-EH1353	Human IFN-β(Interferon Beta) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1356	Human IL-1α(Interleukin 1 Alpha) ELISA Kit	3.91-250pg/mL	2.35 pg/mL
FY-EP1380	Porcine AP17(Apelin 17) ELISA Kit	31.25-2000 pg/mL	9.14 pg/mL
FY-EH1381	Human TNFSF14 ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1387	Human IL-3(Interleukin 3) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1399	Human IL-23(Interleukin 23) ELISA Kit	39.07-2500pg/mL	23.44 pg/mL
FY-EP1420	Porcine BNP(Brain Natriuretic Peptide) ELISA Kit	31.25-2000 pg/mL	9.27 pg/mL
FY-EM1434	Mouse ACTH(Adrenocorticotrophic Hormone) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1436	Human PRL(Prolactin) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL

Cytokine Series Products

Cata NO.	Product Name	Detection Range	Sensitivity
FY-ER1019	Rat ANG2(Angiopoietin 2) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EM1094	Mouse CA9(Carbonic Anhydrase IX) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1106	Human CTSV(Cathepsin V) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1232	Human GP39(Glycoprotein 39, Cartilage) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EH1233	Human CNTF(Ciliary Neurotrophic Factor) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1237	Human F3/TF(Tissue Factor) ELISA Kit	7.82-500pg/mL	4.69 pg/mL
FY-EH1245	Human GCP-2(Granulocyte Chemotactic Protein 2) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1285	Human EPCR(Endothelial Protein C Receptor) ELISA Kit	0.63-40ng/mL	0.38 ng/mL
FY-EH1287	Human FAS/CD95(Factor Related Apoptosis) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-ER1289	Rat FAS/CD95(Factor Related Apoptosis) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-ECH1300	Chicken Neopterin ELISA Kit	0.16-10 ng/mL	0.053ng/mL
FY-EH1341	Human HGF(Hepatocyte Growth Factor) ELISA Kit	125-8000pg/mL	75 pg/mL
FY-EH1396	Human IL-22(Interleukin 22) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM1400	Mouse IL-23(Interleukin 23) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1411	Human EREG(Epregulin) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1417	Human REN(Renin) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH1421	Human CTSC(Cathepsin C) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1447	Human IL-12(Interleukin 12) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1457	Human HLA-B27(Human Leukocyte Antigen B27) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EM1483	Mouse FN(Fibronectin) ELISA Kit	1.57-100ng/mL	0.94 ng/mL
FY-EH1488	Human PC III(Procollagen III) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1490	Human P III NP(N-Terminal Procollagen III Propeptide) ELISA Kit	23.44-1500pg/mL	14.06 pg/mL
FY-EH1494	Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EM1504	Mouse C5a(Complement Component 5a) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EH1512	Human Receptor I for the Fc Region of Immunoglobulin G ELISA Kit	0.79-50ng/mL	0.47 ng/mL
FY-EH1513	Human NCAD(Neural Cadherin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-ER1515	Rat NCAD(Neural Cadherin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-ECH1522	Chicken IL17(Interleukin 17) ELISA Kit	15.63-1000 pg/mL	25.4 pg/mL
FY-EH1523	Human LUM(Lumican) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1525	Human PRR4(Proline Rich Protein 4, Lacrimal) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1526	Human PSAP(Prosaposin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1542	Human TNFRSF1A ELISA Kit	7.82-500pg/mL	4.69 pg/mL
FY-EH1545	Human IL-15(Interleukin 15) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1548	Human 5-LO(Arachidonate 5-Lipoxygenase) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1557	Human ARO(Aromatase) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1573	Human ADAM9(A Disintegrin And Metalloprotease 9) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EH1574	Human ADAM with Thrombospondin Type 1 Motif 4 ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EH1614	Human BMP-3(Bone Morphogenetic Protein 3) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1616	Human FLG-Ab(Filaggrin Antibody) ELISA Kit	1.57-100ng/mL	0.94 ng/mL
FY-EH1619	Human CAPS(Calcyphosine) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1623	Human AFP-L3(Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA	0.63-40ng/mL	0.38 ng/mL
FY-EH1625	Human FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1628	Human CCNB1(Cyclin-B1) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH1629	Human AGC(Aggrecan) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH1631	Human Total Advanced Glycosylation End Product Specific Receptor	78.13-5000pg/mL	46.88 pg/mL

Hot Sale Products

Cata NO.	Product Name	Detection Range	Sensitivity
FY-ER5652	Rat hs-INS(high-sensitivity Insulin) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH5346	Human sFlt-1 ELISA Kit(Replacement)	31.25-2000pg/mL	18.75 pg/mL
FY-EM4950	Mouse NSE(Neuron Specific Enolase) ELISA Kit	78.13-5000pg/mL	46.88 pg/mL
FY-EH3006	Human PGF(Placental Growth Factor) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-ER3164	Rat PI3K2α(Phosphoinositide 3 Kinase 2 Alpha) ELISA Kit	0.16-10ng/mL	0.1 ng/mL
FY-EH3213	Human VF(Visfatin) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EH5272	Human LTF/LF(Lactoferrin) ELISA Kit	0.32-20ng/mL	0.19 ng/mL
FY-EB5834	Bovine Pentraxin-related protein PTX3(PTX3) ELISA kit	0.32-20ng/mL	0.16ng/mL
EY-EB5835	Bovine procalcitonin (PCT) ELISA Kit	31.25-2000pg/mL	18.75pg/mL
FY-EB5831	Bovine Serum amyloid A (SAA) ELISA Kit	6.25-400ng/mL	2.72ng/mL
FY-EB5819	Bovine Tumor Necrosis Factor Alpha(TNF-alpha) ELISA Kit	15.63-1000pg/mL	6.4pg/mL
FY-EB5833	Bovine Interleukin 1-beta(IL1-BETA)ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH5234	Human PLA2G2A(Phospholipase A2, Group II A) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-ER4805	Rat TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-ER2007	Rat cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EU1198	E1(Estrone) ELISA Kit	31.25-2000 pg/mL	8.86 pg/mL
FY-EH5193	Human BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit	31.25-2000ng/mL	18.75 ng/mL
FY-EH4633	Human CRP(C-Reactive Protein) ELISA Kit	0.78-50ng/mL	0.47 ng/mL
FY-EH2008	Human IL-7 ELISA Kit	7.82-500pg/mL	4.69 pg/mL
FY-EH1545	Human IL-15 ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH1559	Human IL-18 ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH3943	Human IL-33 ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH3099	Human IL-34 ELISA Kit	93.75-6000pg/mL	56.25 pg/mL
FY-EH1310	Human FGF21(Fibroblast Growth Factor 21) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EU1013	TGF-β1 ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EH3210	Human VASPIN ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EH1797	Human APLN(Apelin) ELISA Kit	62.5-4000pg/mL	37.5 pg/mL
FY-EH5840	Human Asprosin ELISA Kit	0.16-10ng/mL	0.10 ng/mL
FY-EH5157	Human Cytotoxic T-Lymphocyte Associated Antigen 4 ELISA	0.16-10ng/mL	0.1 ng/mL
FY-EH1405	Human TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	7.81-500pg/mL	3.91 pg/mL
FY-EH4665	Human IL-17A(Interleukin 17A) ELISA Kit	7.81-500pg/mL	4.69 pg/mL
FY-EH1393	Human IL-13(Interleukin 13) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EH4070	Human RNLS(Renalase) ELISA Kit	7.82-500ng/mL	4.69 ng/mL
FY-ER5025	Rat NFκB-p105(Nuclear factor NF-kappa-B p105 subunit)	0.16-10ng/mL	0.1 ng/mL
FY-EH5269	Human hs-CRP(high-sensitivity C-Reactive Protein) ELISA	15.63-1000pg/mL	9.38 pg/mL
FY-EH3015	Human High Sensitive Oxytocin Receptor (OXTR) ELISA	0.16-10ng/mL	0.1 ng/mL
FY-EB5844	Bovine NPY(Neuropeptide Y) ELISA Kit	31.25-2000pg/mL	18.75pg/mL
FY-EM1403	Mouse IFN-γ(Interferon Gamma) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM1406	Mouse TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	7.82-500pg/mL	4.69 pg/mL
FY-EM1445	Mouse IL-1β(Interleukin 1 Beta) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM4663	Mouse IL-2(Interleukin 2) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL
FY-EM1357	Mouse IL-1α(Interleukin 1 Alpha) ELISA Kit	0.79-50pg/mL	0.47 pg/mL
FY-EH1333	Human GDF15(Growth Differentiation Factor 15) ELISA Kit	23.44-1500pg/mL	14.06 pg/mL
FY-EM5315	Mouse IL-6(Interleukin 6) ELISA Kit	31.25-2000pg/mL	18.75 pg/mL
FY-EM4666	Mouse IL-17A(Interleukin 17A) ELISA Kit	15.63-1000pg/mL	9.38 pg/mL